

Product no **AS10 866**

## Horse Ig fraction

### Product information

**Host** | Horse**Purity** | Ig fraction obtained by a proprietary 2-step procedure, in 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2.**Format** | Liquid, 0,2 µm filtered.**Quantity** | 10 mg**Storage** | Store at 2-8 °C.**Additional information** | Horse immunoglobulin fraction is provided in 10 mM sodium phosphate, 0,15 M sodium chloride, pH 7,2, 0,005 % sodium azide is added as preservative

### Application information

||| For applications or usage on species others than stated above Agrisera offers a **payment-after-testing option**. To learn more about this or for any questions on this product, please use the LiveChat option in the left menu bar or contact us at [\*\*email protected\*\*](mailto:support@agrisera.com)