Agrisera Part of Olink® Group \$

This product is for research use only (not for diagnostic or therapeutic use)

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS21 4670 Alcohol dehydrogenase, Saccharomyces cerevisiae

Product information

Immunogen	Alcohol dehydrogenase isolated and purified from UniProt: P00330
Host	Rabbit
Clonality	Polyclonal
Purity	Purified IgG in PBS, pH 7.2.
Format	Lyophilized
Quantity	10 mg
Reconstitution	For reconstitution add 1 ml of sterile, distilled water
Storage	Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to the cap or sides of the tube.

Application information

Expected apparent MW	36,8 kDa (further processed into mature form)
Confirmed reactivity	Saccharomyces cerevisiae
Predicted reactivity	Species of your interest not listed? Contact us

1/1

contact: support@agrisera.com