

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS09 483

AtCCaP1 | vacuolar calcium-binding protein-related

Product information

Immunogen KLH-conjugated synthetic peptide derived from Arabidopsis thaliana CCaP1m UniProt: Q9SXE9, TAIR: At1g62480

Host Rabbit

Clonality Polyclonal

Purity Serum

Format Lyophilized

Quantity 100 μl

Storage

Store lyophilized/reconstituted at -20 °C; make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to the cap or sides of

Additional information

0.1 % sodium azide is added as preservative. For antibody re-suspending information check the tube lable.

Antibodies will detect target protein in a few µg of a crude preparation loaded per well. If purified preparations of vacuolar and plasma membranes are used, one µg load per well should be sufficient.

Application information

Recommended dilution 1:8000 (ELISA), 1:1000 (WB)

Expected | apparent 16.6 kDa

MW

Confirmed reactivity Arabidopsis thaliana

Predicted reactivity Arabidopsis thaliana

Not reactive in No confirmed exceptions from predicted reactivity are currently known

Additional information Protein or membrane sample should be treated at 70°C for 10 min before loading on the gel.

> Diluted antibody solution can be used 2 to 3 times within one month if it contains 0.1 % sodium azide as preservative and is stored at -20°C to -80°C.

Selected references Ide et al. (2007). Transcriptional Induction of Two Genes for CCaPs, Novel Cytosolic Proteins, in Arabidopsis thaliana in the Dark. Plant Cell Physiol. 1:54-65.

Application example



5 μg (1,3) and 2.5 μg (2,5) of lysate from E.coli cells expressing CCaP1 protein/lane were separated on 12 % SDS-PAGE and blotted 1h to PVDF membrane (40 min. at 10 V using BioRad semidry transfer). Filters were blocked 1h with 5 % low-fat milk powder in TBS-T (0.05% Triton X.100). Membranes were washed 5 times with TBS-T, each time in a fresh polystyrene box and probed with anti-CCaP1 antibodies (AS09 483, 1:1000, 1h) and secondary anti-rabbit (1:2000, 1 h). All steps were performed in RT with agitation.