

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS09 601 c-Myc (polyclonal)

Product information

Immunogen KLH-conjugated peptide, sequence MEQKLISEEDLNE, human c-Myc UniProt: Q6LBK7

Host Chicken

Clonality Polyclonal

Purity Immunogen affinity purified IgY in PBS pH 8.5. Contains 0.02% azide.

Format Lyophilized

Quantity 100 μg

Storage Store lyophilized antibodies at -20°C and reconstituted antibodies at 4°C. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to the cap or sides of the tube.

Application information


Recommended dilution 1 : 2000-1 : 5000 (ICC)

Expected | apparent Depends on target protein which is tagged with this epitope

Confirmed reactivity c-Myc

Predicted reactivity c-Myc

Not reactive in No confirmed exceptions from predicted reactivity are currently known



HEK293 cells were transfected with the indicated plasmid (SARS-CoV-2 N myc tagged, PMID: 34799561, TBEV NS3-HA, PMID: 29321318, Rab5 mcherry, Addgene: 4920, GFP-HIS) using genejuice transfection reagent (EMD Millipore) according to the manufacturer's instructions. After 24 hours of transfection, cells were fixed in 4% formaldehyde and permeabilized in PBS containing 0.5% Triton X-100 and 20 mM glycine. Then, cells were stained with the primary anti-Myc tag antibodies at a concentration of 1 μg/mL for 1 hour at room temperature. Followed by three washes in PBS. Cells were then stained using secondary antibodies, goat anti-chicken Alexa488 (1:500, Invitrogen A11039) in PBS containing 2% BSA for 1h at RT. Nuclei were stained using DAPI (1 μg/mL). Images were acquired using a Leica SP8 Laser Scanning Confocal Microscope with a 63x oil objective (Leica) and Leica Application Suit X software (LAS X, Leica).

Courtesy of Dr. Anna K Överby, Molecular Infection Medicine Sweden (MIMS), Section of Virology. Department of Clinical Microbiology Umeå University, Sweden