

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS19 4339 Cyc1 | Cyclin-B1-1

Product information

Immunogen Recombinant Cyclin-B1-1 protein derived from *Arabidopsis thaliana* sequence, amino acids: 1-428. UniProt:

P30183 TAIR: <u>AT4G37490</u>

Host Rabbit

Clonality Polyclonal

Purity >95%, Protein G purified to a total immunoglobulin G fraction.

Format Liquid

Quantity 100 μg

Storage Store at -20 °C or -80 °C, avoid repeated freeze-thaw cycles. Make aliquots to avoid repeated freeze-thaw cycles.

Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material

adhering to the cap or sides of the tube.

Additional information Preservative: 0.03% Proclin 300. Preparation contains: 50% Glycerol, 10 mM PBS, pH 7.4

Reactivity of this antibody on endogenous extract remains to be confirmed.

Application information

Recommended dilution 1:1000-1:5000 (WB)

Expected | apparent 48,4 | 49 kDa

Predicted reactivity | Brassica campestris, Capsella rubella, Eutrema salsugineum, Noccaea caerulescens

Species of your interest not listed? Contact us

Not reactive in No confirmed exceptions from predicted reactivity are currently known

application example



40 and 20 μg of recombinant *Arabidopsis thaliana* CYCB1-1 protein were loaded/well and separated on 10 % SDS-PAGE and blotted 1h to PVDF (0.45 μm) using semi-dry transfer. Blot was blocked with 5 % milk for 2h/RT with agitation. Blot was incubated in the primary antibody at 1.8 μg/ml for 1h/RT with agitation in PBS-T. The antibody solution was decanted and the blot was rinsed briefly twice, then washed 4x in PBS-T at RT with agitation. Blot was incubated in matching secondary antibody (anti-rabbit lgG horse radish peroxidase conjugated) diluted to 1:50 000 in for 1h/RT with agitation. The blot was washed as above and developed with chemiluminescent detection reagent according to manufacture's instructions.