

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS15 3087

Endo-beta-1,4 xylanase (hyperthermophilic enzyme)

Product information

Immunogen His-tagged, recombinant from endo-beta-1,4 xylanase (hyperthermophilic enzyme) of *Thermotoga neapolitana*, overexpressed in *E.coli*, UniProt: <u>B9K775</u>

Host Rabbit

Clonality Polyclonal

Purity Serum

Format Lyophilized

Quantity 50 μl

Reconstitution For reconstitution add 50 μl of sterile water

Storage Store lyophilized/reconstituted at -20 °C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to

the cap or sides of the tube.

Application information

Recommended dilution 1:3000 (WB)

Expected | apparent

MW 119 kDa

Confirmed reactivity Thermotoga neapolitana

Predicted reactivity Thermotoga neapolitana

Not reactive in No confirmed exceptions from predicted reactivity are currently known