

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS10 866

Horse Ig fraction

Product information

Host Horse

Purity | Ig fraction obtained by a proprietary 2-step procedure, in 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2.

Format Liquid, 0.2 µm filtered.

Quantity 10 mg

Storage Store at 2-8 °C.

Additional information Horse immunoglobulin fraction is provided in 10 mM sodium phosphate, 0.15 M sodium chloride, pH 7.2, 0.005 % sodium azide is added as preservative

Application information

||| For applications or usage on species others than stated above Agrisera offers a **payment-after-testing option**. To learn more about this or for any questions on this product, please use the LiveChat option in the left menue bar or contact us at ****