

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS19 4269

PA200 | Proteasome activator PA200

Product information

Immunogen PA200 antibodies were generated against a partial fragment encompassing residues 840 –1140 *Arabidopsis thaliana* protein sequence UniProt: F4JC97-1, TAIR: At3q13330

Host Rabb

Clonality Polyclonal

Purity Serum

Format Lyophilized

Quantity 50 μl

Reconstitution For reconstitution add 50 μl, of sterile water

Storage Store lyophilized/reconstituted at -20°C (short tem, months) or at -80°C (long term, years); once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid

any losses that might occur from material adhering to the cap or sides of the tube.

Application information

Recommended dilution 1:500 (WB)

Expected | apparent

202,9 | ca, 200 kDa

Not reactive in No confirmed exceptions from predicted reactivity are currently known

Additional information Recommended protein load is 20 µg/well, primary antibody dilution: 1: 1000 ON/4°C and optimisation based on

obtained result regarding signal/noise ratio,

Selected references Book et al. (2010). Affinity purification of the Arabidopsis 26 S proteasome reveals a diverse array of plant proteolytic

complexes. J Biol Chem. 2010 Aug 13;285(33):25554-69. doi: 10.1074/jbc.M110.136622.