

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

## Product no AS18 4200-1ml

## Pectic polysaccharide, alpha-1,5-arabinan (monoclonal, clone LM13)

## **Product information**

Immunogen Pectic polysaccharide, alpha-1,5-arabinan, This antibody was isolated from a high throughput screen of many antibodies generated by immunization with a pectic fraction,

**Host** Rat

Clonality Monoclonal

Subclass/isotype IgM

**Purity** Cell culture supernatant.

Format Liquid

Quantity 1 ml

Storage

Store at +4°C (short term) and at -20°C (long term). Make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from any material adhering to the cap or sides of the tube.

**Additional information** 

Contains 0.05% Sodium Azide.

Antibody recognition of arabinans increases with arabinofuranosidase action.

Binds to a specific subset of pectic arabinans, and to longer stretches of 1,5-linked arabinosyl residues that are likely to be more abundant in unbranched arabinans.

## **Application information**

Recommended dilution 1:10 (ELISA, IF)

Confirmed reactivity Higher plants, ferns and mosses

Not reactive in No confirmed exceptions from predicted reactivity are currently known

Selected references

<u>Verhertbruggen</u> et al. (2009). Developmental complexity of arabinan polysaccharides and their processing in plant cell walls. Plant J. 2009 Aug;59(3):413-25.doi: 0.1111/j.1365-313X.2009.03876.x.

Moller et al. (2008). High-throughput screening of monoclonal antibodies against plant cell wall glycans by hierarchical clustering of their carbohydrate microarray. inding profiles Glycoconj J. 2008 Jan;25(1):37-48. doi:

10.1007/s10719-007-9059-7.