

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS03 035-DL650 SPS | Sucrose phosphate synthase, global, DyLight® 650 conjugated (40 μg) Product information

Immunogen	KLH-conjugated synthetic peptide derived from conserved region within plant SPS protein sequences, including <i>Arabidopsis thaliana</i> isoforms 1F Q94BT0, 2F, 3F and 4F. <i>Oryza sativa</i> Q67WN8, <i>Solanum tuberosum</i> Q43845
Host	Rabbit
Clonality	Polyclonal
Purity	Immunogen affinity purified serum, in PBS pH 7.4, conjugated to DyLight® 650.
Format	Liquid
Quantity	40 µg
Storage	Store at 4°C for 12-18 months. A preservative may be added for long time storage, up to 2 years. Shortly spin the tube before use.
	DyLight® 650 has Amax = 652 nm, Emax = 672 nm. DyLight® is a registered trademark of Thermofisher Inc., and its subsidiaries.
Application information	
Recommended dilution	To be determined by end user
Expected apparent MW	120 120-130 kDa (fragments of 30/90 kDa may be detected)
Confirmed reactivity	Arabidopsis thaliana, Colobanthus quitensis Kunt Bartl, Hordeum vulgare, Lycopersicum esculentum, Lycopersicum penelli, Solanum tuberosum, Triticum aestivum, Pinus strobus, Zea mays
Predicted reactivity	Brassica napus, Citrus sinensis, Glycine max, Nicotiana tabacum, Oryza sativa, Physcomitrella patens, Populus balsamifera, Robinia pseudoacaci, Ricinus communis, Saccharum officinarum, Solanum lycopersicum, Theobroma cacao, Vicia faba, Vitis vinifera

Species of your interest not listed? Contact us

- **Not reactive in** No confirmed exceptions from predicted reactivity are currently known.
- Additional information Peptide used to elicit anti-SPS antibodies is perfectly conserved in all isoforms of SPS in plants.
 - Selected references To be added when available, antibody released in May 2023.